

WSI Viewer Quick Guide

Initial viewing, channel display, annotations, and export

1 Open and navigate

- Choose an image in **Sample images**.
- Drag to pan; scroll or pinch to zoom.
- Use viewer controls to reset or enter fullscreen.
- Collapse side panels when more image space is needed.

2 Adjust channels

- Channel controls change presentation only, not source pixels.
- Use checkboxes to show or hide channels.
- Choose LUT color and adjust opacity.
- Use **Reset display** or **Recompute auto** as needed.

3 Work with annotations

- Use **Annotations** and **Names** to show or hide overlays.
- Select one annotation to view or edit its name.
- Enter or leave the name field to commit; Escape cancels a draft.
- Use drawing tools to create annotations and the available delete control to remove one.

4 Protect recent annotation work

- Saving is asynchronous and currently has no visible Saved indicator.
- After creating, moving, renaming, or deleting, pause several seconds before switching images, refreshing, or closing the page.
- If possible, confirm persistence by returning to the image before ending important work.

5 Export

- **Selected annotation** exports the selected region.
- **Entire View** currently means the visible viewer region, not necessarily the whole native slide.
- Native exports are limited to **16,000,000 pixels**. Zoom in and retry if no file is produced.
- A save-location prompt may appear. Confirm the downloaded PNG opens.

6 Environment and safety

- Production normally has no environment warning banner.
- DEVELOPMENT or STAGING banners mean validation only and not clinical use.
- Use only images appropriate for the assigned environment.
- Never copy sensitive production material into development or staging.

Report: blank/distorted images, missing or changed annotations, failed exports, login errors, or unusual delays. Report the image name and approximate time; do not send sensitive image data outside the approved support process.

Before reporting a display problem, try one normal refresh - but never refresh while a recent annotation change may still be saving.